
i-Port Advance™

Injection port • Port d'injection • Injektionsport •
Aplicador de inyecciones • Injctiepoort • Accesso per
iniezioni sottocutanee • Injeksjonsport • Injektioportti •
Injektionsport • Injektionsport • Porta de injeção •
Injekcijski nastavak • Injekční port • Injekciós bemenet •
Port iniekcijny • Enjeksiyon portu • Injekčný port •
Θύρα έγχυσης • Инъекционный порт • 注射端口 •
פתח הזרקה • منفذ الحقن

EN Insertion Sites

FR Sites d'insertion

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FI Asetuskohdat

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DA Indstikssteder

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SL Mesta vstavitve

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SK Miesta zavedenia

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















Indications for Use

The i-Port Advance™ injection port is indicated for patients who administer or receive multiple daily subcutaneous injections of physician prescribed medications, including insulin. The device may remain in place for up to three days to accommodate multiple injections without the discomfort of additional needle sticks. The i-Port Advance may be used on a wide range of patients, including adults and children. The injection port is indicated for single use.

Intended use

The injection port is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The injection port has been tested to be compatible with insulin approved for subcutaneous infusion.

Product Description

The i-Port Advance injection port is a subcutaneous injection port, which combines the injection port with an aid for insertion. The injection port is an integrated part of the inserter. Two injection port models are available: a 6 mm or 9 mm cannula length. The i-Port Advance injection port is a sterile unit, delivered assembled for single use.

Contraindications

- Do not use the same i-Port Advance injection port for longer than three days.
- Do not attempt to reapply the i-Port Advance injection port after use.
- The i-Port Advance injection port is neither intended nor indicated for use with blood or blood products.
- Do not use the i-Port Advance injection port for any purpose other than the indications specified by a healthcare professional.
- Do not use a needle thicker than 28 gauge or thinner than 32 gauge, when injecting into the i-Port Advance. Doing so may damage the septum and may cause incorrect medication delivery.
- Do not use a needle longer than 5/16 in. (8 mm) when injecting into the i-Port Advance injection port. Longer needles may damage the device causing unnecessary punctures of the skin, or tearing or puncturing of the soft cannula, which may result in unpredictable medication delivery.
- Do not use a needle shorter than 3/16 in. (5 mm) when injecting into the i-Port Advance injection port. Shorter needles may not pass through the port septum, which may prevent the medication from properly entering the body.

Warnings and Precautions

- Carefully read through instructions prior to using the i-Port Advance injection port. Failure to follow instructions may result in pain or injury.
- The i-Port Advance injection port is only sterile and non-pyrogenic if the insertion device is undamaged or unopened. Do not use if the package is already open or is damaged, or if the protective liner has fallen off. Ensure the sterility by checking the sterile paper and tamper-proof seal.
- If the i-Port Advance injection port is not securely placed in the inserter with the needle pointing straight ahead prior to insertion, pain or injury may occur.
- Be sure that the needle guard is removed before insertion.
- Do not use if introducer needle is bent or broken.
- Do not inject multiple medications contraindicated for subcutaneous administration into the same insertion site.
- i-Port Advance is a single-use device and should be disposed of immediately after use. Do not clean or re-sterilize.
- Do not reuse the port. Reuse of the port may cause site irritation or infection and may lead to inaccurate medication delivery.
- Incorrect medication absorption, infection, or site irritation may result from improper application to the body or maintenance of the i-Port Advance injection port insertion site, or both. If the insertion site becomes irritated or infected, remove and discard the device and apply a new device to a different location on the body.
- The soft cannula may crimp or bend when applying or while wearing. A crimped or bent cannula may limit or prevent medication delivery through the i-Port Advance injection port. Remove and discard the device immediately if a crimped or bent cannula is suspected. Indications that the i-Port Advance injection port has a crimped or bent cannula include:
 - difficulty pressing down on the plunger of the syringe when injecting
 - injected medication escaping from the device onto the top of the septum
 - injected medication leaking onto the skin or saturating the tape
 - moisture present inside the transparent body of the device
- Make sure that the soft cannula does not extend beyond the introducer needle. If it does, do not use the device. This may result in improper insertion, causing a bent cannula and incorrect medication delivery.
- Do not reinsert the introducer needle into the i-Port Advance injection port. Reinsertion could cause puncturing or tearing of the soft cannula which would result in unpredictable medication delivery or injury.

- Replace the i-Port Advance injection port if the adhesive tape becomes loose or is displaced from its original site. Since the cannula is soft, it may not cause pain if it slips out, and this may occur without notice. The soft cannula must always be completely inserted to receive the full amount of medication from the injection.
- Refer to a healthcare professional for instructions about site rotation. Improper site rotation can result in incorrect medication delivery, infection or site irritation.
- Wash hands, and swab skin with alcohol or similar disinfectant where the patient intends to apply the i-Port Advance injection port. Let skin air dry completely prior to applying the device to the disinfected area.
- Ensure the septum is clean before every use (use purified water).
- Do not put perfumes or deodorants on the i-Port Advance injection port as these may affect the integrity of the device.
- Do not inject into the i-Port Advance injection port through clothing.
- Never point the loaded inserter toward any body part where insertion is not desired.
- Avoid applying the i-Port Advance injection port to sites that contain scar tissue, lipodystrophy, or stretch marks. Avoid sites that are constrained by clothing or accessories, or that are subjected to rigorous movement during exercise. Avoid sites that are under a belt or on the waistline.
- Do not inject more than 75 times through a single device.
- The i-Port Advance injection port may have up to 2.60 μL or 0.26 units of medication that remain in the port. This can delay or prevent your medication from entering the body. Consult with the healthcare professional of your prescribed medication to determine if this impacts therapy management.

Recommendations


- Choose insertion sites as recommended by a healthcare professional. The choice of site depends on treatment and patient-specific factors such as body composition and physical activity level. Refer to a healthcare professional for instructions about insertion site rotation.
- Test blood glucose per a healthcare professional's instructions to ensure that insulin delivery at the insertion site is appropriate.
- Test blood glucose 1 to 3 hours after using the i-Port Advance injection port to ensure insulin delivery performs as anticipated, and measure blood glucose on a regular basis as specified by a healthcare professional.
- Avoid applying and using a new i-Port Advance injection port prior to bedtime unless blood glucose can be checked 1 to 3 hours after injection.

- In the case of an unexplained high blood glucose reading, remove and discard the i-Port Advance injection port and apply a new device to a different location on the body since the soft cannula could be dislodged or partially clogged. Continue to monitor blood glucose readings to ensure insulin is absorbed. Should suggested actions not provide a resolution to the problem or should blood glucose remain high, contact a physician or healthcare professional.
- Contact a healthcare professional if symptoms do not match test results and all of the instructions in this manual have been followed.
- Consult with a healthcare professional if there are any questions or concerns regarding diabetes management.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to a national authority.

Storage and disposal

- Store injection ports in a cool, dry location at room temperature. Do not store injection ports in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid the risk of needle stick.
- Dispose of the used injection port according to local regulations for biohazardous waste.

i-Port Advance™ injection port

- | | |
|------------------------|------------------------|
| (A) Sterile liner | (H) Lid |
| (B) Center of inserter | (I) Round indentations |
| (C) Inserter | (J) Septum |
| (D) Finger grips | (K) Adhesive tape |
| (E) Introducer needle | (L) Protective liner |
| (F) Injection port | (M) Soft cannula |
| (G) Needle guard | |

Insertion

Unfold the front cover of this booklet to view images that are numbered to the steps below.

- 1** Wash hands with soap and water.
Select an insertion site on a recommended area (shown in gray).
- 2** Clean the insertion site with a disinfectant. Be sure the area is air dry before inserting i-Port Advance injection port. If needed, remove hair around the insertion site to prevent improper attachment of the adhesive tape to the skin. Contact a healthcare professional if there are any questions regarding proper insertion site cleaning.
- 3** Pull the red tab and remove the seal.
- 4** Remove the sterile liner.
- 5** Securely hold the inserter and remove the lid.
- 6** Remove the protective liner from the tape.

- 7** Gently twist to remove the needle guard from the introducer needle. Check that the soft cannula does not exceed the introducer needle.
- 8a** Place fingers on the finger grips and hold the inserter so that the introducer needle is facing down or away from the body.
- 8b** With the other hand, pull up on the center of the inserter until it is locked in an upright position.
- 9** Place the inserter against the skin and place fingers on the round indentations. Squeeze both round indentations at the same time to insert the injection port.
- 10** Press down on the center of the inserter to secure the injection port to the skin.
- 11** Grip the center of the inserter and gently pull to remove the inserter from the injection port.
- 12** Gently massage the adhesive tape onto the skin.
- 13** Put the lid back on the insertion device and dispose of the inserter according to local requirements.

Using the i-Port Advance injection port

- 1** Clean the septum with purified water before each use.
- 2** Puncture the septum using an injection needle.
- 3** Slowly inject medication into the injection port.
- 4** Remove the injection needle.

Removal

- 1** Carefully loosen the edges of the adhesive tape, and remove the injection port slowly by pulling straight out. Make sure that the needle is still attached to the adhesive and pulled completely out of the subcutaneous tissue (skin). If the needle is missing, immediately contact a healthcare professional for assistance.

© 2020 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (TM) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries. i-Port Advance™

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Medtronic Österreich GmbH
Medtronic B.V.