Infusion set • Equipo de infusión • Perfuseur • Infusionsset • Set di infusione • Infusieset • Infusjonssett • Infusionsset • Zestaw infuzyjny • Infuzní souprava
• Conjunto de infusão • ערכת העירוי • Σετ έγχυσης • Infuzijos rinkinys • 주입 세트 • 注射器 • مجموعة التسريب • Инфузионный набор

dana inset™ II

Instructions for use • Instrucciones de uso • Mode d'emploi • Anwendungshinweise • Istruzioni per l'uso • Gebruiksaanwijzing • Bruksanvisning
• Bruksanvisning • Instrukcja użycia • Návod k použití • Indicações de utilização • הוראות השימוש • Οδηγίες χρήσης • Naudojimo instrukcija
• 사용 설명서 • 使用说明 • إرشادات الاستعمال • Инструкции по применению

SOOIL

- en** Use-by (Expiry date) YYYY-MM-DD
- es** No utilizar después de (fecha de caducidad) (AAAA-MM-DD)
- fr** Date limite d'utilisation (date d'expiration) AAAA-MM-JJ
- de** Verwendbar bis (Verfallsdatum) JJJJ-MM-TT
- it** Utilizzare entro (data di scadenza) AAAA-MM-GG
- nl** Gebruiken tot (vervaldatum) JJJJ-MM-DD
- no** Brukes innen (utløpsdato) ÅÅÅÅ-MM-DD
- sv** Utgångsdatum ÅÅÅÅ-MM-DD
- pl** Zużyć przed (data ważności) RRRR-MM-DD
- cs** Použitelnost (expirační doba) RRRR-MM-DD
- pt** Utilizar até (Data de validade) AAAA-MM-DD
- he** תאריך תפוגה (לשימוש עד) DD-MM-YYYY
- el** Χρήση έως (ημερομηνία λήξης) EEEE-MM-HH
- lt** Naudoti iki (galiojimo pabaigos data) MMMM-MM-DD
- ko** 사용기한(만료일자) YYYY-MM-DD
- zh** 使用期限(有效期) 年-月-日
- ar** تاريخ انتهاء الصلاحية: يوم-شهر-سنة
- ru** Срок использования (срок действия) ГГГГ-ММ-ДД

- en** Do not re-use
- es** No reutilizar
- fr** Ne pas réutiliser
- de** Nicht wiederverwenden
- it** Non riutilizzare
- nl** Niet hergebruiken
- no** Må ikke brukes om igjen
- sv** Får ej återanvändas
- pl** Nie używać ponownie
- cs** Nepoužívejte opakovaně
- pt** Não reutilizar
- he** לא לשימוש חוזר
- el** Μην επαναχρησιμοποιείτε
- lt** Nenaudoti pakartotinai
- ko** 재사용하지 마십시오
- zh** 请勿重复使用
- ar** لا تعيد استخدامه
- ru** Для одноразового использования

- en** Sterilized using ethylene oxide
- es** Esterilizado con óxido de etileno
- fr** Stérilisé à l'oxyde d'éthylène
- de** Mit Ethylenoxid sterilisiert
- it** Sterilizzato con ossido di etilene
- nl** Gesteriliseerd met ethyleenoxid
- no** Sterilisert med etylenoksid
- sv** Steriliserad med etylenoxid
- pl** Sterylizowany przy użyciu tlenku etylenu
- cs** Sterilizováno etylenoxidem
- pt** Esterilizado com óxido de etileno
- he** ממוסטר מראש במצעות אתילן-אוקסיד
- el** Αποστειρώθηκε με χρήση αιθυλενοξειδίου
- lt** Sterilizuota naudojant etileno oksidą
- ko** 산화 에틸렌을 이용한 멸균처리
- zh** 使用环氧乙烷消毒
- ar** تم التعقيم باستخدام أكسيد الإيثيلين
- ru** Стерилизовано оксидом этилена

- en** Caution: see instructions for use
- es** Precaución: ver instrucciones de uso
- fr** Attention : voir notice d'utilisation
- de** Achtung: siehe Gebrauchsanweisung
- it** Attenzione: consultare le istruzioni per l'uso
- nl** Let op: zie de gebruiksaanwijzing
- no** OBS! Se bruksanvisningen
- sv** Försiktighet! Mer information finns i bruksanvisningen
- pl** Przestroga: patrz instrukcja użytkowania
- cs** Upozornění: viz návod k použití
- pt** Atenção: consultar as instruções de utilização
- he** הרתוח: יש ליינין בהוראות השימוש
- el** Προσοχή: ανατρέξτε στις οδηγίες χρήσης
- lt** Dėmesio: žr. naudojimo instrukciją
- ko** 주의: 사용 설명서 참조
- zh** 注意事项: 见使用说明
- ar** تحذير: راجع إرشادات الاستعمال
- ru** Внимание: см. инструкцию по применению

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Indications for use:

The Dana inset™ II infusion set is indicated for the subcutaneous infusion of insulin administered by an external pump. The infusion set is indicated for single use.

INTENDED USE

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

DESCRIPTION

Dana inset™ II is an infusion device which combines an infusion set with an aid for insertion (serter). As the infusion set is an integrated part of the serter, Dana inset™ II is delivered assembled and ready for use as a single, sterile unit.

CONTRAINDICATIONS

Dana inset™ II is neither intended nor indicated for use in intravenous infusion of blood or blood products.

WARNINGS

- Carefully examine the instructions for use before inserting Dana inset™ II. Failure to follow instructions may result in pain or injury.
- Dana inset™ II is only sterile and non-pyrogenic if the insertion device is undamaged and unopened. Do not use if the insertion device is already open or has been damaged. Ensure sterility by checking if the sterile paper and tamper-proof seal are undamaged.
- When using Dana inset™ II for the first time, do so in the presence of a healthcare provider. Since different individuals have varying amounts of subcutaneous layer, ensure that an appropriate cannula length is chosen to fit your needs.
- Inaccurate medication delivery, site irritation or infection may result from improper insertion or maintenance of your infusion site.
- Do not re-insert the introducer needle into the infusion set. This could cause a tearing or damage of the soft cannula and may lead to unpredictable medication flow.
- Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.
- Be sure that the needle guard is removed

before insertion is made.

- If the soft cannula bends during insertion, apply a new Dana inset™ II immediately at a new site.

PRECAUTIONS

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- When unwinding, release the tubing with caution as a hard pull of the tubing can result in damage to the infusion device. Ensure that the infusion set is properly in place in the serter when the tubing is fully released – see picture 8.
- Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult pump manufacturer's instructions for use.
- Dana inset™ II is a single use device and should be disposed of immediately after use. Do not clean or re-sterilize.
- Do not in any way put disinfectants, perfumes, deodorants, cosmetics or other products containing disinfectants in contact with the connector or the tubing. This may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose or is displaced from its original site. Since the cannula is soft, it will not cause

any pain if it becomes dislodged, and this may take place without notice. The soft cannula must always be completely inserted to receive the full amount of medication.

- Replace the infusion set every two to three days, or per your healthcare provider's instructions.
- If your insertion site becomes irritated or infected, replace the infusion set and use a new site until the first site has healed.
- Wash your hands before disconnecting and reconnecting the infusion set. Consult your healthcare provider on how to compensate for missed medication when disconnected.
- Never point a loaded insertion device towards any body part, where insertion is not desired.

RECOMMENDATIONS

- Check your blood glucose level one to three hours after insertion of the Dana inset™ II. Measure your blood glucose on a regular basis. Discuss this with your healthcare provider.
- Do not change your infusion set just prior to bedtime, unless blood glucose can be checked one to three hours after insertion.
- If your blood glucose level becomes unexplainably high or an occlusion alarm

occurs, check for clogs and leaks. If in doubt, change your infusion set since the soft cannula could be dislodged or partially clogged. Discuss plan to bring down your blood glucose level with your healthcare provider should any of these problems arise. Test your blood glucose level to ensure problem is corrected.

- Carefully monitor your blood glucose levels when disconnected and after reconnecting.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.



STORAGE & DISPOSAL

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid risk of needle stick.
- Dispose the used infusion set according to local regulations for biohazardous waste.

Dana inset™ II - Preparation and insertion

Read all instructions carefully before using the Dana inset™ II infusion set. For information about pump therapy, possible errors, and potential risks and their causes, consult your pump manufacturer's instructions for use before connecting Dana inset™ II. When using Dana inset™ II, follow proper hygiene procedures.

Please see the corresponding illustrations on page 80-82.

CONTENTS

- | | |
|---------------------|-----------------------|
| [A] Lid | [F] Connector needle |
| [B] Needle guard | [G] Tubing |
| [C] Soft cannula | [H] Introducer needle |
| [D] Adhesive tape | [I] Insertion device |
| [E] Cannula housing | [J] Sterile paper |

IMPORTANT NOTE!

- [K] Three raised dots – press here when lifting the lid off (see Figure 6)
- [L] Lined indentations – press here when loading Dana inset™ II (see Figure 10)
- [M] Round Indentations – press here simultaneously when releasing Dana inset™ II (see Figure 12)
- [N] Disconnect cover (see Figure B)

- 1 Wash your hands before inserting Dana inset™ II.
- 2 Suggested insertion sites. Choose an insertion site as recommended by your healthcare provider. Do not use an area right next to your previous infusion site.
- 3 Clean the area for insertion with a disinfectant as directed by your healthcare provider. Be sure the area is air dry before inserting Dana inset™ II.
- 4 Pull the red tape to remove the seal.
- 5 Remove the sterile paper.
- 6 Press the three raised dots on each side of the lid firmly with one hand and lift up the lid with the other hand.

CAUTION: Be careful not to bend or touch the Dana inset™ II introducer needle when preparing the set and infusion set.

IMPORTANT NOTE! If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead prior to insertion, pain or minor injury may occur. To avoid this, please ensure that the following steps are taken gently so the infusion set is not loosened / removed by accident.

- 7 **Unwinding the tubing:** Carefully twist the beginning of the tubing out of the slot. Now

unwind the tubing from the insertion device by pulling it carefully upwards.

CAUTION: Do not touch the Dana inset™ II introducer needle when unwinding the tubing. Be careful not to pull too hard when the tubing is fully released. You could accidentally remove the infusion set from the introducer needle.

- 8** Connect Dana inset™ II to a pump reservoir and the pump. Prime the infusion set. Ensure there are no air bubbles in the tubing. Consult pump manufacturer's instructions for use.

CAUTION: When priming the infusion set, hold Dana inset™ II with the needle pointing down to ensure that medication does not come in contact with the adhesive backing paper.

- 9** Pull up gently to remove the adhesive backing paper.
- 10 Prepare the inserter:** Place fingers on lined indentations. Press lined indentations on each side. Pull spring until you hear a "click".
- 11** Carefully remove the needle guard by gently twisting and then pulling it. Check that the soft cannula does not exceed the introducer needle.

CAUTION: Before positioning Dana inset™ II over a clean insertion site, place the tubing in the slot to ensure that the tubing is not caught

underneath the infusion set during insertion.

- 12 Insert Dana inset™ II:** Position Dana inset™ II over the insertion site. Simultaneously press the round indentations on either side of the insertion device to insert Dana inset™ II.

- 13** Push gently on the center of the inserter to secure the adhesive on the skin.

- 14** Remove the inserter/introducer needle by grasping the center and pulling it gently back. Gently massage the adhesive tape onto the skin.

Fill the empty space in the cannula with the following amount of insulin:

U-100 insulin:

6 mm : 0.1 units (0.001 ml)

9 mm : 0.2 units (0.002 ml)

- 15** Push the lid back in place until you hear a "click".

- A** Dana inset™ II allows you to temporarily disconnect your pump without having to change the infusion set. Hold the adhesive in place by placing a finger just in front of the cannula housing, and gently squeeze the sides of the connector. Pull the connector needle straight out from the cannula housing.

- B** Insert the disconnect cover into the cannula housing until you hear a "click".

- C** If necessary, prime the infusion set until medication dispenses from the needle.

CAUTION: When priming the infusion set, hold Dana inset™ II with the needle pointing down.

- D** Place a finger on the adhesive tape just in front of the cannula housing.

- E** Push the connector needle straight in until you hear a "click".

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SOOIL Development Co., Ltd.