Infusion set ▪ Equipo de infusión ▪ Set de perfusion ▪ Infusionsset ▪ Set di infusione ▪
İnfüzyon seti ▪ Infusieset ▪ Infusionssæt ▪ مجموعة التسريب ▪ Infusionssett ▪ Infusionsset ▪
Σετ έγχυσης ▪ Infuzní souprava ▪ Infuusionsetti ▪ Zestaw infuzyjny ▪
Set za infuziju ▪ Conjunto de infusão ▪ סט ניירוי ▪ Infūzijas komplekts ▪
Система за инфузия ▪ Infusioonikomplekt ▪ Infúziós szerelék ▪ Innrennslisett ▪
Infuzinis rinkinys ▪ Set de perfuzie ▪ Инфузионный набор ▪ Infúzna súprava ▪ Infuzijski set

Instructions for use ▪ Instrucciones de uso ▪ Mode d'emploi ▪ Gebrauchsanweisung ▪ Istruzioni per l'uso ▪
Kullanım talimatları ▪ Gebruiksaanwijzing ▪ Brugsanvisning ▪ إرشادات الاستخدام ▪ Bruksanvisning ▪ Bruksanvisning ▪
Οδηγίες χρήσης ▪ Návod k použití ▪ Käyttöohje ▪ Instrukcja użycia ▪ Upute za korištenje ▪ Instruções de utilização ▪
הוראות השימוש ▪ Lietošanas instrukcija ▪ Инструкции за употреба ▪ Kasutusjuhend ▪
Használati útmutató ▪ notkunarleiðbeiningar ▪ Naudojimo instrukcijos ▪ Instrucțiuni de utilizare ▪
Инструкции по применению ▪ Návod na použitie ▪ Navodila za uporabo

en
es
fr
de
it
tr
nl
da
ar
no
sv
el
cs
fi
pl
hr
pt
he
lv
bg
et
hu

See page 10
Véase página 14
Voir page 18
Siehe Seite 22
Vedere pagina 26
Bkz. sayfa 30
Zie pagina 34
Se side 38
راجع الصفحة 42
Se side 46
Se sida 50
Βλ. σελ. 54
Viz str. 58
Katso sivu 62
Patrz strona 66
Pogledajte stranicu 70
Ver página 74
ראה בעמוד 78
Skatiet 82. lpp.
Вижте стр. 86
Vt lk 90
Id. 94. oldal

is
lt
ro
ru
sk
sl

Sjá bls. 98
Žr. 102 psl.
A se vedea pagina 106
См. стр. 110
Pozrite si stranu 114
Glejte stran 118

REF

| | |
|----|-------------------------------------|
| en | Reference / model number |
| es | Referencia/número de modelo |
| fr | Référence/numéro de modèle |
| de | Bestellnummer/Modellnummer |
| it | Riferimento/Numero di modello |
| tr | Referans/model numarası |
| nl | Referentie- / modelnummer |
| da | Reference/modelnummer |
| ar | الرقم المرجعي / رقم الموديل |
| no | Referanse-/modellnummer |
| sv | Referens-/modellnummer |
| el | Αναφορά / αριθμός μοντέλου |
| cs | Referenční číslo / číslo modelu |
| fi | Viite-/mallinnumero |
| pl | Numer katalogowy/numer modelu |
| hr | Kataloški broj / broj modela |
| pt | Referência/número de modelo |
| he | מספר דגם / מסמך |
| lv | Atsauces/modela numurs |
| bg | Референтен номер/номер на модел |
| et | Viite-/modelinumber |
| hu | Hivatkozási szám / típuszám |
| is | Tilvisunarnúmer/gerðarnúmer |
| lt | Nuoroda / modelo numeris |
| ro | Număr de referință/model |
| ru | Код / номер модели |
| sk | Referenčné číslo/číslo modelu |
| sl | Referenčna številka/številka modela |

| | |
|----|--|
| en | Caution: see instructions for use |
| es | Precaución: ver instrucciones de uso |
| fr | Attention : voir le mode d'emploi |
| de | Achtung: Gebrauchsanleitung durchsehen |
| it | Attenzione: consultare le istruzioni per l'uso |
| tr | Dikkat: Kullanma talimatlarına bakın |
| nl | Let op: zie de gebruiksaanwijzing |
| da | Forsigtig: Se brugsanvisningen |
| ar | تحذير: راجع إرشادات الاستخدام |
| no | OBS: se bruksanvisningen |
| sv | Försiktighet! Mer information finns i bruksanvisningen |
| el | Προσοχή: ανατρέξτε στις οδηγίες χρήσης |
| cs | Upozornění: viz návod k použití |
| fi | Vaara: katso käyttöohjeet |
| pl | Przeostroga: patrz instrukcja użycia |
| hr | Oprez: pogledajte upute za korištenje |
| pt | Atenção: consultar as instruções de utilização |
| he | היזהרו: ראו הוראות השימוש |
| lv | Uzmanību! Skatiet lietošanas instrukciju |
| bg | Внимание: вижте инструкциите за употреба |
| et | Hoiatus! Lugege kasutusjuhendit |
| hu | Vigyázat: olvassa el a használati útmutatót |
| is | Varúð: sjá leiðbeiningar um notkun |
| lt | Atsargiai: žr. naudojimo instrukciją |
| ro | Atenție: vedeți instrucțiunile de utilizare |
| ru | Внимание: см. инструкцию по применению |
| sk | Upozornenie: pozrite si návod na použitie |
| sl | Pozor: glejte navodila za uporabo |

| | |
|----|----------------|
| en | Manufacturer |
| es | Fabricante |
| fr | Fabricant |
| de | Hersteller |
| it | Produttore |
| tr | Üretici |
| nl | Fabrikant |
| da | Producent |
| ar | الشركة المصنعة |
| no | Produsent |
| sv | Tillverkare |
| el | Κατασκευαστής |
| cs | Výrobce |
| fi | Valmistaja |
| pl | Producent |
| hr | Proizvođač |
| pt | Fabricante |
| he | יצרן |
| lv | Ražotājs |
| bg | Производител |
| et | Tootja |
| hu | Gyártó |
| is | Framleiðandi |
| lt | Gamintojas |
| ro | Producător |
| ru | Изготовитель |
| sk | Výrobca |
| sl | Proizvajalec |

| | |
|----|----------------------------------|
| en | Date of Manufacture YYYY-MM-DD |
| es | Fecha de fabricación AAAA-MM-DD |
| fr | Date de fabrication AAAA-MM-JJ |
| de | Herstellungsdatum JJJJ-MM-TT |
| it | Data di produzione AAAA-MM-GG |
| tr | Üretim Tarihi YYYY-AA-GG |
| nl | Datum van fabricage JJJJ-MM-DD |
| da | Produktionsdato ÅÅÅÅ-MM-DD |
| ar | تاريخ التصنيع: يوم-شهر-سنة |
| no | Produksjonsdato ÅÅÅÅ-MM-DD |
| sv | Tillverkningsdatum ÅÅÅÅ-MM-DD |
| el | Ημερομηνία κατασκευής ΕΕΕΕ-MM-HH |
| cs | Datum výroby RRRR-MM-DD |
| fi | Valmistuspäivämäärä VVVV-KK-PP |
| pl | Data produkcji RRRR-MM-DD |
| hr | Datum proizvodnje GGGG-MM-DD |
| pt | Data de fabrico AAAA-MM-DD |
| he | תאריך הייצור DD-MM-YYYY |
| lv | Izgatavošanas datums GGGG-MM-DD |
| bg | Дата на производство ГГГГ-ММ-ДД |
| et | Tootmiskuupäev AAAA-KK-PP |
| hu | Gyártás ideje (ÉÉÉÉ-HH-NN) |
| is | Framleiðsludagur ÁÁÁÁ-MM-DD |
| lt | Pagaminimo data MMMM-MM-DD |
| ro | Data de fabricație AAAA-LL-ZZ |
| ru | Дата изготовления ГГГГ-ММ-ДД |
| sk | Dátum výroby RRRR-MM-DD |
| sl | Datum proizvodnje LLLL-MM-DD |

| | |
|----|---|
| en | Global Trade Item Number |
| es | Número de artículo de comercio mundial |
| fr | Code article international (GTIN) |
| de | Globale Artikelnummer |
| it | Codice GTIN (Global Trade Item Number) |
| tr | Global Ticari Öğe Numarası |
| nl | Globale artikelnummer |
| da | Globalt handelsvarenummer |
| ar | رقم الصنف التجاري العالمي |
| no | Globalt handelsartikkelnummer |
| sv | Globalt handelsnummer |
| el | Αριθμός παγκόσμιου εμπορικού στοιχείου |
| cs | Globální číslo obchodní položky |
| fi | Kansainvälinen kauppatavaranumero |
| pl | Globalny numer jednostki handlowej |
| hr | Globalni broj trgovačke jedinice |
| pt | Número de item de comércio global |
| he | מספר פריט סחר גלובלי |
| lv | Globālās tirdzniecības preces numurs |
| bg | Глобален номер на търговската единица |
| et | Globaalne kaubaartikli number |
| hu | Globális cikkszám |
| is | Alþjóðlegt vörunúmer |
| lt | Visuotinės prekybos prekės numeris |
| ro | Numărul global al articolului comercial |
| ru | Глобальный номер единицы товара |
| sk | Globálne číslo obchodnej položky |
| sl | Globalna trgovinska številka izdelka |

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Indications for use:

The product is indicated for subcutaneous infusion of medication administered by an external pump.

INTENDED USE

The product is intended for the subcutaneous infusion of immunoglobulins for the treatment of primary immune deficiency, apomorphine and foslevodopa/foscarbidopa for parkinson's disease, morphine (hydromorphone, morphine sulfate and morphine chloride) for pain management therapy, saline for mild to moderate dehydration and deferoxamine mesylate for thalassemia. The product has been tested to be compatible with the above-mentioned drugs approved for subcutaneous infusion. Clinical evaluation supports infusion set use from age 2 years and up.

PRODUCT DESCRIPTION

The product (drawing 1a, 1b, 1c, page 12) consists of a cannula housing (H) with a 90-degree soft cannula (I) and a tubing (J) with a standard pump connector (luer-lock) (L) combined with an integrated insertion device (A). The product is delivered ready for use as a sterile unit.

CONTRAINDICATIONS

The product is neither intended nor indicated for intravenous (I.V.) infusion of medication, including blood and blood products.

CLINICAL BENEFITS

- Continuous subcutaneous infusion to mimic stimulation of receptors (Parkinson), serum IG range (Immunodeficiency), iron chelation treatment (Thalassemia), prevent or treat mild to moderate dehydration, and pain relief (Pain Management).
- Less invasive administration route associated with fewer adverse events and more locations to infuse compared to the intravascular route.
- Less injections (insertions) compared to injection therapy.

WARNINGS

- **Carefully examine the instructions for use before inserting the product. Failure to follow instructions may result in pain or injury.**
- **Do not use the product if the box and/or blister pack are already opened or damaged. Make sure that the safeguard (B) on top of the insertion device (A) has not been removed. If this occurs please use a new product.**
- **Do not reposition the product as this may compromise the adhesion of the adhesive tape (E).**
- **If blood is seen in the cannula housing (H) and/or tubing (J), change to a new product and a new site.**
- **If the adhesive tape (E) becomes loose, change to a new product and a new site.**
- **The product integrity may be compromised if exposed to solvents containing substances,**

e.g. disinfectants, perfumes and deodorants and/or cosmetics and skincare products.

- **Do not reuse the infusion set. Reuse of the infusion set may cause infection, site irritation, or damage to the cannula/needle. A damaged cannula/needle may lead to inaccurate medication delivery.**
- **Never attempt to repair a clogged tubing (J) while the product is connected to the body. An uncontrolled quantity of medication may accidentally be infused. If clogged, change to a new product and a new site.**
- **Choice of product and insertion site is closely linked and should always be done by a healthcare professional.**
- **Never point the insertion device towards any body part where insertion is not desired.**

PRECAUTIONS

- When using the product for the first time, do so in the presence of a healthcare professional.
- The product is a single use item which must be disposed of after use. Do not clean or re-sterilize.
- Wash hands with soap and water when temporarily disconnecting the product. Seal the cannula housing (H) and site connector (K) with the covers (B and M) provided. Consult with a healthcare professional on how to compensate for missed medication when disconnected, and for how long the product can be disconnected from the pump.
- Do not leave air in the tubing (J) when priming.

Make sure the medication is kept and handled according to the instructions for the medication.

- Avoid contact with used soft cannula (I).
- Do not change to a new or unknown product without consulting a healthcare professional for correct handling.
- Change the product in consultation with a healthcare professional following the below recommendations:
 - after 12 hours of use with apomorphine,
 - after 12 hours of use with immunoglobulins,
 - after 18 hours of use with deferoxamine mesylate,
 - after 72 hours of use with foslevodopa/foscarbidopa,
 - after 72 hours of use with morphine,
 - after 72 hours of use with saline.
- If using a new insertion device and cannula housing (1b), follow illustrations 1 to 10 in this instruction. Prime and connect the tubing (J) separately referring to illustrations G, H and I.

STORAGE & DISPOSAL

- Store infusion sets in a cool, dry location at room temperature (15-25°C). Do not store infusion sets in direct sunlight, or in high humidity.
- Dispose of the insertion device in a proper sharps container according to local laws.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

ADDITIONAL INFORMATION FOR THE PATIENT

- Use an insertion site as recommended by your healthcare professional. Do not use an insertion site which is infected, swollen, red, irritated, bruised or is bleeding. Do not insert in an area with lumps, hard tissue or bumps. Be careful not to insert into bone or muscle as this may cause pain, damage the product and affect medication absorption. Choose an area that does not have constrictive clothing or other objects that can pull on the infusion set.
- Please refer to your healthcare professional for insertion site rotation and general recommendations.
- Check the insertion site frequently to ensure that the soft cannula (I) remains firmly in place. If it is not in place, change to a new product and a new site. The soft cannula (I) must always be completely inserted to ensure delivery of the full

amount of medication. If the insertion site becomes red or inflamed, change to a new product and a new site.

- The product is available with various soft cannula (I) and tubing (J) lengths. Follow the recommendations of your healthcare professional regarding appropriate soft cannula (I) and tubing (J) length.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.


ADDITIONAL INFORMATION FOR THE HEALTHCARE PROFESSIONAL

- Please consider treatment and patient specific factors when choosing an insertion site.
- The choice of soft cannula length should be based on treatment and patient specific factors such as physiology and activity level. A too long soft cannula may result in insertion pain and/or insertion into a bone or muscle and a too short soft cannula in leakage or irritation at the insertion site. In both cases the infusion must be terminated and treatment should not be resumed until an infusion set with another soft cannula length is inserted at a new insertion site.
- Please inform the patient about the correct choice of insertion site in relation to the choice of soft cannula length.


Preparation and insertion

Go through the instructions for use carefully before inserting the product, as failure to follow instructions may result in inaccurate medication delivery, infection and/or site irritation. For information about pump therapy, including connection and priming procedures, possible errors, and potential risks and their causes, consult your pump manufacturer's instructions for use before using the product. When using the product for the first time, do so in the presence of a healthcare professional.


THE PRODUCT (1a)

INSERTION DEVICE & CANNULA HOUSING (1b)

TUBING (1c)

1 CONTENTS

[A] Insertion device

[B] Safeguard (Optional: transparent cover for cannula housing)

[C] Activation button (for insertion of the infusion set)

[D] Indication for tubing direction

[E] Adhesive tape

[F] Paper

[G] Insertion site window

[H] Cannula housing

[I] Soft cannula

[J] Tubing

[K] Site Connector

[L] Pump connector (luer-lock)

[M] White cover for site connector

- Wash your hands with soap and water.
- Use an insertion site as recommended by your healthcare professional. The choice of site depends on treatment and patient specific factors. Check that the insertion site is free of skin irritation such as redness, scar tissue, etc.
- Clean the area of insertion with a disinfectant as recommended by your healthcare professional. Be sure the area is air dried before inserting the product. If needed, remove hair around the insertion site to ensure proper attachment of the adhesive tape (E) to the skin.
- Open the blister pack and take out the contents.
- Pull gently to remove the paper (F) from the adhesive tape (E). Be careful not to touch the adhesive tape (E).
- Remove the safeguard (B) from the insertion device (A) by gently squeezing the sides of the safeguard (B) while pulling it straight out.
- Before placing the insertion device (A) against the skin, note the arrows (D) on the insertion device (A) showing where the tubing (J) will be connected after insertion.
- Gently stretch the skin and ensure that the surface is flat and taut. Position the insertion device (A) on the skin and press the red activation button (C) completely down to insert the infusion set. The insertion is complete when the insertion device (A)

is released from the cannula housing (H) with adhesive tape (E). The product is now in safe mode with the activation button (C) pressed down. The product is now inserted.

- 10** Press the adhesive tape (E) onto the skin. Change to a new product and a new site, if the adhesive tape (E) becomes loose.
- 11** Remove the paper from the tubing (J).
- 12** Connect the tubing (J) to the pump or syringe driver by turning the pump connector (L) clockwise.
- 13** Remove the white cover (M) from the site connector (K).
- 14** Prime* the tubing (J) according to the instructions provided by the pump manufacturer. Be sure there is no air in the tubing (J).
- 15** Place a finger on the cannula housing (H) while pushing the site connector (K) straight in until you hear a “click”.
- 16** Frequently check the insertion site to ensure that the soft cannula (I) remains in place and that the site is not irritated. Change to a new product and a new site if necessary. Remember to store the covers (B and M) for later use. Dispose of trash according to local recommendations.

DISCONNECTING THE TUBING

- A** Wash your hands with soap and water.
- B** Place a finger on the cannula housing (H) and gently squeeze the sides of the site connector (K). Pull the site connector (K) straight out from the cannula housing (H).
- C** Place the white cover (M) on the site connector (K).
- D** Insert the transparent cover (B) into the cannula housing (H) until you hear a “click.”

RECONNECTING THE TUBING

- E** Wash your hands with soap and water.
- F** Remove the transparent cover (B) by placing a finger on the cannula housing (H) and gently squeezing the sides of the cover. Pull the cover straight out of the cannula housing (H).
- G** Remove the white cover (M) from the site connector (K).
- H** Make sure there is medication flow and no air in the tubing (J).
- I** Place a finger on the cannula housing (H) while pushing the site connector (K) straight in until you hear a “click”.

REMOVAL OF THE INFUSION SET

- J** Wash your hands with soap and water.
- K** Place a finger on the cannula housing (H) and gently squeeze the sides of the site connector (K). Pull the site connector (K) straight out from the cannula housing (H).
- L** Remove the infusion set by carefully loosening the adhesive tape (E) around the cannula housing (H) before pulling the soft cannula (I) out of the skin.
- M** Disconnect the tubing (J) from the pump or syringe driver by turning the pump connector (L) counter clockwise. Dispose of the infusion set according to local recommendations. For available sharps containers, please visit your local pharmacy.

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