



The Esteem+ portfolio

The Esteem+ portfolio offers a range of one-piece convex products with various depths to help provide comfort, freedom, security and reassurance to your patient, all designed for different stoma shapes and sizes.



Esteem+ Convexity Variant Guide*



Full assessment of the stoma and abdomen should be made when selecting convexity variant**


Esteem+ SOFT

Firm/toned abdomen

Key benefits:

- Shallow wafer depth minimises pressure on delicate skin
- Designed to apply minimal pressure and help conform to the body
- Suitable for use immediately post-operatively

Suited to protruding and flush stomas




Esteem+ FLEX

Skin folds & creases

Key benefits:

- Fits to your body contours for comfort and easy positioning
- Hydrocolloid skin barrier to keep the pouch firmly in place

Suited to protruding, flush and retracted stomas




Esteem+ FIRM

Soft /flaccid abdomen

Key benefits:

- Good seal for patients with soft flaccid abdomen
- Encourages protrusion of flush or retracted stomas

Suited to flush and retracted stomas



*Based on Convexity Plateau Variants Guidelines / Considerations selection guide.

**Patients should be assessed by their HCP/ Stoma Nurse for appropriate product choice and use of convexity.

Data on file Convatec Inc.

To find out more about our ostomy solutions, please contact our Customer Care Team

Australia: Freecall 1800 335 276 or visit www.convatec.com.au

New Zealand: Call 0800 225 4309 or visit www.convatec.com/en-nz/

ALWAYS FOLLOW THE DIRECTIONS FOR USE.

® / TM indicate trademarks of the Convatec group of companies

© 2024 Convatec Inc. AP-70107-AUS-ENG-v1 (v1.4) O568.3 September 2024



convatec
— forever caring —